

From: Steve Ross 2-9510, skross @ anl.gov  
To: Detector Pool  
Date: 12/07 onward...  
Subject: Sarnoff CAM1M100 Camera

## Sarnoff CAM1M100 High Frame Rate CCD camera.

### **REVISION HISTORY:**

Initial version 2/14/2008

This memo provides the specifications for the APS Detector Pool high frame rate visible light Sarnoff CCD camera. It also is a User's Guide for operation.



### **1.0 ADVANTAGES OF THIS CAMERA:**

A charge coupled device (CCD) camera relies on transferring charge from an imaging pixel to the output amplifier by a so-called bucket brigade, moving the photo-charge one row of pixels to the next (vertical transfer), then shifting a given row out through a readout amplifier (horizontal transfer). This readout amplifier can have higher and higher bandwidth, but at the expense of electronic noise. Hence after a certain point, a better

way to achieve fast frame readout is to increase the number of output amplifiers, while keeping each at a lower readout speed for lower noise. The Sarnoff series is one of the main commercial cameras to take this approach, having 16 output amplifier ports. Hence when you look at an image, you will see 16 sub-regions in the raw data. It is never possible to make all regions exactly the same light-to-signal gain, exactly the same dc offset. It is straight forward to correct this in later software.

**Compared with other visible light CCD cameras, the Sarnoff camera has a very good quantum efficiency (thus raising sensitivity), low noise, true 12 bit dynamic range at 10 millisecond readout. It is a full-frame CCD, the entire area is part of the imaging area (not interline, not frame transfer.) Hence it typically requires a shutter.**

## 2.0 DETECTOR POOL TAG NUMBERS:

DP00244	Sarnoff Fast visible CCD camera
DP00245	Sarnoff fast shutter
DP00246	Sarnoff fast shutter controller

## 3.0 SPECIFICATIONS:

**Format:** Silicon Charge coupled device (CCD) 1024 x 1024 pixels, back side thinned.

**Optical Lens Mount:** C mount

**Pixel Size:** 16 um x 16um

**Visible light quantum efficiency:** about 70%, see curve below,

There is no window in front of the CCD in the camera housing

**Number of output ports:** 16

**Readout time:** 10 milliseconds full frame.

**Typical Operating temperature:** Cooled by 2 stage thermo electric coolers to -30C

**Electronic Noise** at 100 frames/sec (essentially no integration time) 20eRMS (see data below).

**Dynamic Range, Linearity, and full well capacity:**

The amount of charge a CCD can hold in its imaging pixel is set by the bias on the CCD, and by the doping of the semiconductor itself. This is an upper bound on the “full well” capacity, but not universally useful because linearity (digital signal out versus photons in) degrades before this limit is reached. The Sarnoff quoted camera pixel full well capacity of 140,000 electrons is instead set by the full scale of the 12 bit analog to digital converter (ADC). That is, the gain on the pixel amplifier chains is such that 140,000 electrons equal 4095 DN. The camera is linear better than <0.26% through the full scale. (see data below).

Note we estimate the dynamic range of CCD (as opposed to the ADC)

$$: \quad \frac{140,000}{25eRMS} = 5600 \text{ dynamic range of CCD}$$

$$12 \text{ bits} = 4096$$

Measured at frame rate is 100 fps full frame, or 200 fps binned by 2x2

**Offset:** 128 DN is the dark offset value as set in the standard camera.

Dark current <1na/cm<sup>2</sup> at 20C, which means <200e/pix/sec at -30C

So with exposures longer than about 15 sec the dark current would be a factor. (See appendix section of dark current noise, if you wish to estimate this for yourself.)

Other definitions: “Back side thinned”. The “front side” of a CCD is defined at that part of the silicon which has the gate structure, the lines on the silicon which brings the clock phases to the pixel. This structure can attenuate light, blue (400 nm) more than red (800 nm). Hence we can thin the back of the silicon wafer down from hundreds of microns to very thin few microns, and call this the photo active side. The advantage is better visible light quantum efficiency, but the expense is a more fragile wafer, one that probably cannot be bonded to anything else mechanically. The Sarnoff CCD is this back side thinned.


## SPECIFICATIONS (TYPICAL)

Resolution	1024 x 1024 pixels
Pixel Size	16 um x 16 um
Aperture	16.4 mm x 16.4 mm
Max. Frame Rate <sup>1</sup>	100 fps or 200 fps binned by 2
Binning	2 x 2 or Programmable
Scanning	Split Full Frame
Camera Outputs	2 x 50 MHz LVDS via MDR26
Data Format	12 bit (4096 DN)
Dynamic Range	3200:1
Antibloom	2:1
Lens Mount	F-mount, C-mount
Camera Size <sup>2</sup> (HWD)	17.2 cm x 11.2 cm x 17.7 cm
Mass <sup>2</sup>	< 1.5 kg
Operating Temp.	10 to 35° C
Power Supply	90-260 VAC, 47-63 Hz
Power Supply Size	16.4 cm x 8.5 cm x 5.5 cm, 0.6 kg
Power Dissipation	< 50 Watts

<sup>1</sup> Thermal management is required at the highest frame rates • <sup>2</sup> Without lens

## 4.0 MEASURED PERFORMANCE DATA:

4.1 Visible light Quantum efficiency curve is given below. Typical phosphors and scintillator emit at about 550 nm, but this does vary.



**4.2 Pixel Defect spec:** No bad columns. Pixel responsivity: Non uniformity in bright frames will be within (at 390 nm Typical) [-15%, +20%]. This applies to all pixels, excluding the first 3 and last 3 columns and first 3 and last 3 lines. The following exceptions are allowed: Up to 500 pixels could be defined as blemishes. A blemish is defined as a pixel, which gives a signal within -50% of the average signal. Any pixel outside the [- 50%,+20%] range will be considered a defect pixel. Up to 50 defect pixels are allowed that fall outside the blemish requirements. For each blemish pixel, the average of 8 Nearest Neighbors will not be less than 45% of the overall average.

The camera consists of four internal printed circuit boards into a fifth mother board, or focal plane array board.


### 4.3 TYPICAL CAMERA NOISE DATA:

This data was taken by the manufacturer, running at the full frame rate (about 100 fps). We provide this as backup to the specifications listed above.

PhotonTransfer

Port 1 - PTC DATA					
Exp Time (ms)	Signal Mean	Signal Noise	Read Noise	Shot Noise	Signal Variance
10	1352.051	6.751	0.763	6.708	45.579
15	1490.445	7.039	0.861	6.998	49.554
20	1634.471	7.443	0.87	7.404	55.399
25	1778.735	7.739	0.733	7.701	59.889
30	1923.315	7.995	0.793	7.958	63.913
35	2067.83	8.176	0.709	8.14	66.844
40	2212.095	8.567	0.787	8.533	73.387
45	2355.018	8.598	0.787	8.564	73.922
50	2499.323	9.06	0.912	9.027	82.079
55	2642.859	9.271	0.823	9.239	85.945
60	2786.001	9.487	0.82	9.456	89.997
65	2930.738	9.672	0.804	9.642	93.551
70	3074.484	9.843	0.803	9.813	96.876
75	3219.383	10.142	0.894	10.113	102.857
80	3364.492	10.402	0.89	10.374	108.206
85	3508.592	10.699	0.783	10.672	114.474
90	3653.248	10.84	0.783	10.813	117.496
95	3797.131	10.933	0.787	10.906	119.52
100	3942.83	10.965	0.735	10.939	120.238


Adjust Best Fit Line for Slope = 1/2	
e-/DN =	23 <= Adjust this
Slope =	0.5 <= Until this is 0.5
Camera System Inputs:	
ASIC Gain =	2
Signal Attenuation =	1
Dark Offset (DN)	5
Gain Constant	
K (e-/DN):	23
Sensitivity	
CCD Sensitivity (uV/e-)	10.8
Read Noise	
Ave Read Noise (DN)	0.8
Ave Read Noise (e-)	19
Dynamic Range (SNR)	
SNR	5042
SNR (dB)	74



**4.4 TYPICAL CAMERA LINEARITY DATA: This data was provided by the manufacturer, and serves as backup to the specifications listed above.**

Non-Linearity	
Non-linearity (%)	0.26

Port 2 - Linearity Data			
LED Pulse (ms)	Signal Mean	Best Fit Line	Deviation
0	1.743	2.49	-0.75
6.66668	197.775	200.73	-2.95
13.33332	396.939	398.97	-2.03
20	596.211	597.21	-0.99
26.66668	794.371	795.44	-1.07
33.33332	993.69	993.68	0.01
40	1192.049	1191.92	0.13
46.66668	1391.217	1390.16	1.06
53.33332	1590.399	1588.40	2.00
60	1789.416	1786.63	2.78
66.66668	1988.993	1984.87	2.12
73.33332	2188.931	2183.11	2.82
80	2384.466	2381.35	3.12
86.66668	2581.472	2579.59	1.89
93.33332	2779.257	2777.82	1.43
100	2976.897	2976.06	0.83
106.66668	3173.06	3174.30	-1.24
113.33332	3369.612	3372.54	-2.93
120	3564.547	3570.78	-6.23
m			
29.73571799			
b			
2.49086775			



## 5.0 USERS GUIDE:

### 5.1 OVERALL CAMERA SYSTEM -

Camera (no lens attached, just the aluminum mounting bracket).




Rear view of camera.




Items of interest:

- (1) Note the trigger IN connector (SMB I think, anyway don't lose this cable!). This goes back to the control/interlock box.
- (2) The DB9 connector on top right is for the TEC, which really can only be used if dry gas is set up on the CCD/lens, and the water chiller.
- (3) There is a connector for the dry gas N2 in (top/middle), and vented out.
- (4) There are water in/out connections on the bottom middle.
- (5) And finally on the bottom there are the two camera link data output ports. Cable these to the frame grabber, wire 1 to input 1, wire 2 to input 2, as labeled.

Sarnoff Power Supply




Water Chiller (on cart)



Interlock system



### E-BRICK



### X-ray Shutter





### Shutter Controller




## APPENDIX. SUPPLEMENTAL INFORMATION.

**THERMO ELECTRIC COOLING SYSTEM.** There is a DB9 connector on the side of the camera which allows the user to input current for the TEC system, and to read out the camera's temperature. See TEC Interface Drawing.pdf, the interface to the thermo electric coolers:

(email 5/26/2005)

- [TEC Interface \(see attached drawing\)](#)
- [TEC data:](#)
- [Mode: spectroscopic \(128\\_bin\)](#)
- [Frame period: 200us](#)
- [Twater: 20C](#)
- [CCDTEMP=1.25V \(0C\); Vtec=1.3V, Itec=0.6A](#)
- [CCDTEMP=0.94V \(-10C\); Vtec=2.6V, Itec=1.2A](#)
- [CCDTEMP=0.62V \(-20C\); Vtec=4.4V, Itec=2.2A](#)
- [CCDTEMP=0.33V \(-30C\); Vtec=7V, Itec=3.6A](#)



**ANL SUPPLIED COMPUTER:**

**SOFTWARE.**

We are currently using Sarnoff supplied software which can only take and store one image at a time. This software needs to be replaced by Brian Tieman-like software. This is under development.

**SOFTWARE AND CONTROL, and FIRMWARE.** Sarnoff will provide a C++ based SDK (**optional**) for software control of the camera system. The camera will have multiple clocking “modes” of operation, including a spectroscopy thin strip exposure mode.

**HARDWARE.**

(Ref: Tad’s email of 6/2/2005). We purchase a computer with faster SCSI disks. It will run Windows XP, but could run LINUX etc. Inside this computer is a camera-link based Coreco frame grabber.

Item / description	Part no.	Unit price	Quantity	Ext.price
» HP Workstation xw4200	Base	\$4,380.01	1	\$4,380.01
Microsoft® Windows®				
HP Workstation xw4200	DU936AV			
Microsoft® Windows® XP Professional	DY527AV#ABA			
Intel® Pentium® 4 3.80GHz/1MB 800FSB w/EM64T	DY536AV			
ATI FireGL V3100 128MB PCI-E	PE950AV			
2GB (4x512) DDR2-533 ECC unb	PB340AV			
146GB U320 SCSI 15K (1st)	DY556AV			
146GB U320 SCSI 15K (2nd)	DY568AV			
16X DVD+/-RW, DL, LightScribe, Win (1st)	PR509AV			
U320 SCSI Controller - no RAID, no external connector	DZ541AV			
Floppy Disk Drive	DY569AV			
HP PS/2 Standard Keyboard	DY577AV#ABA			
HP USB Optical Scroll Mouse	DY581AV			
HP xw4200 Localization kit	DY526AV#ABA			
Symantec Norton Antivirus 2004	DY588AV#ABA			
<b>Warranty</b>				
3/3/3 (parts/labor/next business day on-site) limited warranty				
<b>Cart Total:</b>		<b>\$4,380.01</b>		

## CORECO FRAME GRABBER (PCI bus, into above computer). X64-CL-Full\_Dsheet\_web\_March04.pdf

### Specifications<sup>1</sup>

<b>Board</b>	Half length PCI 3.1 64-bit 66MHz compliant 5V and 3.3V slot compatible
<b>Acquisition</b>	Acquisition rates up to 660MB/s Horizontal Size (mb/s/max): 6 bytes/255MB Vertical Size (mb/s/max): Line scan cameras: 1 line to infinity Area scan cameras: 1 line to 16 million/frame variable length frames Onboard frame buffer memory up to 32B (32MB standard) One 1024 x 1024 or 1024 x 255 input lookup table <sup>2</sup> Single slot resolution supports 1 Medium, 1 Base or 2 synchronized Base Camera Link cameras Interfaces to digital area scan or line scan color or monochrome cameras Supports standard multiTap Camera Link configurations: 1 Full camera 8 x 8-bit 1 Medium camera: 4 x 8-bit, 4 x 10-bit, 4 x 12-bit, 1 x 20-bit/RGB and 1 x 26-bit 1 Base camera: 2 x 8-bit, 2 x 10-bit, 2 x 12-bit, 1 x 14-bit, 1 x 16-bit, and 1 x 24-bit/RGB Alternate multiTap Camera Link configurations support: 4 x 14-bit, 4 x 16-bit, 1 x 48-bit or 1 x 64-bit <b>Pixel Formats</b> Monochrome 8, 10, 12, 14, 16, or 26-bit/RGB <b>Transfers</b> Realtime transfers to system memory: PCI-32 bus: 32 bits @ 33MHz PCI-64 bus: 64 bits @ 66MHz PCI-X bus: 64 bits @ 66MHz On-the-fly tap adjustments for multiple tap area scan and line scan cameras
<b>Controls</b>	Comprehensive event notification includes: start/end-of-frame, sequence or P/line events One independent TTL/VDCS trigger input programmable as active-high or low (edge or level trigger) One strobe TTL output for area scan and line scan cameras One PC independent RS-232 COM port provides seamless interface to MS Windows applications Quadrature (AQ) shaft-encoder inputs for external web synchronization; supports cascaded divide and multiply Up to 12 interrupt driven I/Os <sup>3</sup> permit external event synchronization
<b>Power Output</b>	Power-on-reset fused +1.2V/+5V DC output at 1.5A
<b>Software</b>	Microsoft Windows NT 4.0, Windows 2000, and Windows XP compliant Full support of the Sapeva programming package Compatible with Microsoft Visual Studio 6.0 and .Net (unmanaged code only), Visual Basic 6.0 and Borland C++ Builder 5.5 or higher
<b>System Requirements</b>	PCI-64 or PCI-32 compliant system and 64MB system memory
<b>Dimensions</b>	5.2" (21cm) Length x 4.50" (11.7cm) Height
<b>Temperature</b>	0° C (32° F) to 55° C (131° F) Relative Humidity: up to 95% (non-condensing)
<b>Markings</b>	PCB class B - approved CE class B - approved

<sup>1</sup> Last updated April 2004

#### Notes:

1. Please contact Coreco Imaging for more information on Pixel Formats
2. For 1 x 8-bit/tap or 1 x 16-bit/tap maximum
3. Optional module requires auxiliary slot

**CORECO**  
**IMAGING**

# SHUTTER

Strictly speaking, the shutter really is the beamline’s “problem”, we do not adequately support this part of CCD cameras yet. To allow testing in my lab, I purchased a Uniblitz x-ray shutter capable of fast bursts. We know that this shutter is not perfect for beamlines, and we can discuss.

<b>XRS6 SHUTTER SPECIFICATIONS</b>	Page 1 of 2	<b>UNIBLITZ®</b> BY VINCENT ASSOCIATES
Vincent Associates 803 Linden Avenue, Rochester, NY 14625 web www.uniblitz.com e-mail vincentassociates@uniblitz.com	toll-free 800.828.6972	

**FEATURES**

- 6mm diameter aperture
- Pt-Ir shutter blade, capable of blocking x-ray energy up to 30KeV.
- Exposure repetition rate continuously variable from DC-50Hz.
- Electronic synchronization system included.
- Activated by an electronic pulse through UNIBLITZ® patented shutter drive systems.
- Non-resonant design allows instantaneous changes in repetition rate and duty cycle.
- No optical surface when open provides 100% transmittance

The UNIBLITZ® XRS6 series is especially designed for x-ray applications. The innovative blade design allows beam extinction >10<sup>4</sup> up to 30KeV x-ray energy. The XRS6 shutter is well suited for specific x-ray applications such as x-ray crystallography and can open within 3.2msec at a maximum rate of 50 exposures per second. Precision control and reliability can be expected as with all UNIBLITZ products.

This shutter is programmable, and is activated by an electronic pulse generated by Vincent’s patented UNIBLITZ drive systems. The shutter will follow this pulse, which allows the user to program the exposure duration and frequency.

The XRS6 has a 6mm clear aperture and equipped with the electronic synchronization system. A vacuum compatible version is also available, please consult customer service for further information.

**ELECTRICAL**


Coil Resistance	12 ohms
Pulse Voltage to Open	+65VDC
Hold Voltage <sup>1</sup>	+5VDC

**MECHANICAL**

Wgt. Un-Cased	1.20 oz (.03 kg)
Wgt. Cased	7.46 oz (.21 kg)
Operating Temp.	0°C to +80°C
Max. Opening Bounce	15%
Max. Closing Bounce	5%
Max. Frequency of Operation (CONT/BURST) <sup>2</sup>	10 Hz / 50 Hz
Number of Blades	1

**TIMING**

Typical timing values (msec.) using UNIBLITZ drive equipment and measured with UNIBLITZ shutters equipped with .010" PtIr shutter blades.



(Timing in msec.)

Q-A Delay time on opening after current is applied	1.2
A-C Transfer time on opening	3.2
O-C Total opening time	4.4
B-F Min. equivalent exp. time	5.6
C-E Min. dwell time with min. input pulse	2.0
E-G Transfer time on closing	3.8
A-G Total window time	5.0

MET: Min. exposure time 6.4  
TEP: Typical exposure pulse >6.4

The question regarding enhancement of shutter speed with the application of user supplied lubricants has been repeatedly asked. It is our experience that lubricating the shutter blades will actually slow the shutter down and eventually render the shutter inoperable. UNDER NO CIRCUMSTANCES SHOULD ANY TYPE OF LUBRICANT BE APPLIED TO THE SHUTTER BLADE AREA.

**PRODUCT OPTIONS**

<b>XRS6S 2 P 0 -100</b>				
APERTURE SIZE	HOUSING	BLADE FINISH	ELECTRONIC SYNCHRONIZATION	MOUNTING OPTIONS
XRS6S - 6mm	1 - UNCASED 2 - #2 CASE	P - .010" THICK PT (10% IRIUM 90% PLATINUM)	0 - OMIT SYNC. 1 - ELECTRONIC SYNC.	21 ZEISS APOVERT TYPE 22 OLD STYLE NIKON TYPE 23 OLYMPUS TYPE 24 OLYMPUS TYPE 25 LEICA TYPE 26 NIKON TYPE 28 OLYMPUS II TYPE 29 NIKON TYPE 29 LEICA TRANSMITTED TYPE 27 NIKON TYPE 27 NIKONCFODAL TYPE 22 NIKON BB TYPE 100 MOUNTING RING 100 C-MOUNT ADAPTER (MALE) 100 C-MOUNT ADAPTER (FEMALE) 110 T-MOUNT ADAPTER

<sup>1</sup>Voltage level required across actuator coil when being held in the open position.  
<sup>2</sup>CONTinuous frequency rating specified at shutter's minimum exposure pulse. BURST frequency rating specified for (4) four seconds maximum with (1) one minute minimum between bursts. Frequency measurements are taken in free air, 25° C ambient, actuator coil equipped with heat sink. For additional information on maximum sustained frequencies obtainable, please contact one of our technical represent-

## **DIGITAL FRAME GRABBER:**

- (1) Coreco X64-CL-FULL, frame grabber, with on-board 2GB, 85 MHZ (see below)
- (2) Coreco Sopera Lite Software
- (3) High end PC Computer, running WindowsXP, with fast SCSI disks, see spec's below.
- (4) Water Chiller, Model Thermo-Electron/Neslab Merlin series M-25, part number 262 112 030 000, high temperature audible alarm, temp range +5 to +35C, cooling 810 watts at 20 C, pump PD-1, 1.5 gpm at 60 psi, temp stability 0.1C, power from 110 VAC.  
Pump: [www.thermo.com/com/cda/product/detail/1,1055,1000001155234,00.html](http://www.thermo.com/com/cda/product/detail/1,1055,1000001155234,00.html)  
Thermo Electron Corp/ Neslab phone (603)430 2297, fax (603)436-8411, Attn: Scott.
- (5) DC power source (greater than 2 amps) for the TE cooler.
- (6) To monitor the temperature Argonne will supply a digital meter (similar to a DVM).

## **OPTICS.**

Need to discuss typical lens tube – microscope objective then imaging lens, essentially a microscope. C mount. Same as CoolSnap, but someday might modify to include a fast shutter. I have researched shutters, and rejected 2 so far, settling on the UniBlitz.

Emil Trakhtenberg modified the tube lens to allow focus. I have the drawings. Tube lens was supplied by Desna Laboratories (914) 948 9502, see “Tube lens quote.PDF”