

Diagnostics Beamline Development

(Sector 35)

Bingxin Yang and Alex H. Lumpkin

The missions of the APS diagnostics beamlines are:

- (1) To support user operations by providing on-line information about the electron beam on a 24-hour basis;
- (2) To support machine physics studies by providing photon diagnostics on beam properties, transverse and longitudinal;
- (3) To conduct research and development of diagnostics for future light sources.

Strategies:

Utilize and develop time-resolved imaging with optical and x-ray synchrotron radiation to support the mission.

ADVANCED PHOTON SOURCE

OUTLINE

I. Overview of Sector 35

II. Beam Quality Assurance for User Operations

- (A) Pinhole camera upgrade: beam position, size, and emittance data tracking
- (B) Bunch purity and bunch cleaning

III. Undulator-Based Electron Beam Measurements


- (A) Emittance measurements
- (B) Relative beam energy measurement
- (C) Energy-dependent e-beam property measurements
- (D) Observation of lowest effective emittance in APS-SR straight sections: 3.5 nm rad

IV. Summary and Future Developments

ADVANCED PHOTON SOURCE

I. Overview of Sector 35

Sector 35 Layout



Three beamlines:

- 35-BM: X-ray pinhole camera, time-resolved imaging (visible light)
- 35-ID: Undulator-based diagnostics
- 36-AM: Pinhole camera at low-dispersion source, bunch purity monitor

ADVANCED PHOTON SOURCE

SECTOR 35 IMPORTANT EVENTS

- Visible light CCD camera 3/18/95
Image of the first APS SR beam
- X-ray pinhole camera (BM: high dispersion point) 3/96
Pinhole in air, detector in tunnel, $M = 0.44$
- X-ray pinhole camera (BM) 4/97
Pinhole in vacuum, detector in tunnel, $M = 1.06$
- ID beam divergence measurement 5/2/98
Monochromator in air, distance @ 28.4 m, 36.5 m
- X-ray pinhole camera (AM: low dispersion point) 6/98
Pinhole in vacuum, detector in hutch, $M = 1.55$
- X-ray streak camera first image (AM) 11/97
- BM pinhole camera improvement 5/98-1/99
Upgrade in-tunnel detector, stabilize water temperature,
improved video processing
- Second video frame grabber implemented 8/98
Mobile MV-200 used at Sector 35
- ID horizontal beam size measurement 8/99
Pinhole slits after monochromator, $M = 0.277$
- ID beam divergence measurement 9/98
Monochromator in vacuum @28.4m
- Simultaneous tracking of beam size and divergence 1/99
35-ID beamline
- ID-based e-beam measurements 1999
Beam energy spread, energy-dependent properties,
low emittance to 3.5 nm·rad

ADVANCED PHOTON SOURCE


II. Beam Quality Assurance for User Operations

(A) Beam position, size and emittance: pinhole camera upgrade and data tracking

The pinhole camera at 35-BM continued to be our workhorse in 1998-1999. Upgrades performed during 1998 - 1999:

- Stabilizing cooling water temperature to $\pm 0.1^\circ\text{F}$ (ASD/MEC)
- Installing new imaging detector (better optics, 80 μm thick scintillator)
- Improving video signal processing

Data quality has been significantly improved. At 1-Hz bandwidth, we have better than 1- μm stability in measured beam size, as shown by the following data set.




The rms horizontal beam size during the current month recorded by the APS data logger: Significant variations can be seen during machine studies on September 1, 7, 14-15, and 21.

ADVANCED PHOTON SOURCE

Insertion device gap changes have strong effect

During user runs, the beam size can vary up to $2.5\ \mu\text{m}$, largely due to ID gap changes by the user. The following data (2-week and 2-day tracking) show strong (anti)correlation of total energy loss by the APS insertion devices with emittance.




ADVANCED PHOTON SOURCE

The beam position, size, and emittance data are available to the Operations Group through EPICS (1 update / sec.)



and to the users through the Web




ADVANCED PHOTON SOURCE


Diagnose subtle beam instabilities

Having a stable measurement of beam size helped identify a subtle transverse beam instability (> 85 mA) and a longitudinal beam instability (> 101 mA) during user runs in 1998. The instabilities were suppressed. (Lumpkin et al. PAC99)


35-BM source size
(horizontal)



35-ID divergence
(horizontal)



35-ID source size
(horizontal)



ADVANCED PHOTON SOURCE

Resolution measurement of pinhole camera


The measured vertical beam size was assumed to be the quadrature sum of the camera resolution and the true beam size

$$\sigma_{\text{exp}}^2 = \sigma_R^2 + \sigma_y^2.$$

For a single bunch fill, the beam lifetime is dominated by Touschek scattering (A. Ropert, LS-275), which is proportional to the bunch volume. Fitting the measured vertical beam sizes and lifetime to the simple function gives:

Pinhole camera resolution = 22 μm

Vertical beam size = 16 μm (smallest observed)



0.2% vertical coupling is the lowest observed to date.

ADVANCED PHOTON SOURCE

Summary for pinhole camera upgrade beam position, size and emittance data tracking

- In FY1999, the performance of the 35-BM x-ray pinhole camera has reached a submicron level of stability in transverse beam size measurements due to hardware and software upgrades.
- The deduced emittance data is now available to the users through the Web (ASD/OPS).
- During user operations, the APS storage ring horizontal beam size at the BM is normally stable within a 2- μm range (1.4% @ BM, ~3% @ ID). Most variations of horizontal beam size and emittance are induced by insertion device gap changes.
- In at least two incidences, the x-ray pinhole camera has enabled the identification and correction of subtle beam instabilities before users realized their existence.
- The camera resolution has reached 22 μm in mid-1999. It was used to directly observe a vertical coupling of 0.2%, the lowest on the APS storage ring so far.

Further upgrades


(needed for operations with 0.1% coupling, $\sigma_y = 10 \mu\text{m}$)

- The extension of the pinhole camera beamline is in progress. It is expected to improve resolution to the level of 10 - 15 μm .
- Approaches for further improvement in resolution to the level of 5 - 10 μm are being studied.
- Improvement of the stability of transverse size measurement to < 0.5 μm level (1 Hz bandwidth) will be pursued in FY2000.

ADVANCED PHOTON SOURCE

(B) Bunch purity and bunch cleaning

A bunch-purity monitor was implemented by cloning the system in Sector 3 (E. Alp)




Schematic of Sector 35 bunch purity monitor

(The APD is located in 35-ID-A station, taking x-rays from 36-AM pinhole line)

ADVANCED PHOTON SOURCE

During user operation, bunch-purity is routinely checked after each fill




Control screen for the bunch-purity monitor


ADVANCED PHOTON SOURCE

Support bunch-cleaning experiments

The bunch-purity monitor was also used to support bunch-cleaning experiments, an effort towards eventual injections without the Positron Accumulator Ring.



Charge distribution before bunch cleaning



Charge distribution after bunch cleaning


ADVANCED PHOTON SOURCE

III. Undulator-based Electron Beam Measurements

(A) Emittance measurement with undulators

Introduction (a simplified picture)


1. Divergence measurement with monochromatic undulator beam



Monochromator select photon energy: $\omega < \omega_1$

$$\text{Beam Divergence} = w / L$$


2. Pinhole camera for source size measurement



$$\text{Source size} = w * S / S'$$

ADVANCED PHOTON SOURCE

Twiss parameters and systematic errors (geometrical optics approximation)



Electron Beam waist location: $Z = \alpha / \gamma$

- Apparent divergence is larger due to the contribution of e-beam size
- Apparent source size is smaller due to the directional nature of undulator radiation

ADVANCED PHOTON SOURCE

Twiss parameters and systematic errors (geometrical optics approximation)

Effective divergence from x-ray beam cross section ($Z = \alpha / \beta$)

$$S_m^2 \sigma_{x',eff}^2 = (S_m - Z)^2 \sigma_x^2 + S_m^2 \sigma_m^2 + \sigma_{x0}^2$$

Effective beam size from pinhole camera (or any small aperture optic)

$$\sigma_{x,eff}^2 = \frac{S^2 \left(\sigma_{x0}^2 + \sigma_\gamma^2 + \frac{Z^2}{\beta_0^2} \sigma_\gamma^2 \right)}{\sigma_{x0}^2 (S - Z)^2 + \sigma_\gamma^2 S^2 + \sigma_{x0}^2} \sigma_{x0}^2$$

Product of the two effective measurements

$$\varepsilon_{x,eff} = \sigma_{x,eff} \sigma_{x',eff} = \frac{\varepsilon}{R}$$

with correction factor

$$\frac{1}{R} = \frac{S}{S_m} \sqrt{\frac{\beta_0^2 + (S_m - Z)^2 + S_m^2 \frac{\sigma_m^2}{\sigma_{x'}^2}}{\beta_0^2 + (S - Z)^2 + S^2 \frac{\sigma_\gamma^2}{\sigma_x^2}}} \cdot \left(1 + \frac{\beta \sigma_\gamma^2}{\beta_0 \sigma_{x'}^2} \right),$$

a slow function of Twiss parameters. The systematic effects are nearly completely cancelled for $S \approx S_m$.

ADVANCED PHOTON SOURCE


More complications: Fresnel diffraction point spread function of pinhole camera

A pinhole camera point spread function (PSF) is best modeled by broadband Fresnel diffraction pattern.

$$f_5(\tilde{x}) = \frac{1}{2\pi\sqrt{N}} \int_0^2 \frac{dt}{\eta t} \left\{ e^{-\frac{\tilde{\varphi}_0^2 \tilde{\sigma}_0^2 + \alpha^2 (\psi_0 + \tilde{\varphi}_0 \tilde{x})^2}{2\eta^2}} \sin \frac{\psi_0 + \tilde{\varphi}_0 \tilde{x}}{\eta^2} + (\tilde{x} \rightarrow -\tilde{x}) \right\}$$


The observed pinhole data display strong Fresnel diffraction signatures.

INTEGRATED VERTICAL INTENSITY PROFILE SHOWS FRESNEL DIFFRACTION SIGNATURE



ADVANCED PHOTON SOURCE

EXPERIMENTAL SET UP



Emittance correction factor, $R \approx 0.92$

The **APS Diagnostics Undulator** is optimized for beam diagnostics

- Short period: $1.8 \text{ cm} \times 198 \text{ period}$
- Low power : 25.9 keV @ 30 mm gap, $< 1 \text{ W}$ x-ray power
- Low divergence : $2.6 \text{ } \mu\text{rad}$
- Low field error: rms slope error $< 1 \text{ } \mu\text{rad}$

ADVANCED PHOTON SOURCE

Monochromator optimized for divergence measurement

- Direct imaging of unmasked undulator beam: Interpretation straightforward
- Single crystal monochromator: Avoid artifacts associated with double crystal monos
- Short distance (0.5 m) from crystal to scintillator screen: Minimize effect of heat-induced crystal distortion
- Thin silicon crystal (0.3 mm) can be used in Bragg (400) or Laue (220) reflection: Allow substantial hard x-ray flux to pass for source size measurement. In vacuum chamber with 1-mm Be window, located @ 27 m from ID
- Detector: YAG crystal 0.5 mm thick + 7:5 optics + CCD camera, located @ 27.5 m

Pinhole camera parameters ($M = 0.32$)

- Pinhole slits: water cooled 6-mm tungsten blades, located @ 28.5 m, $25\ \mu\text{m} \times 25\ \mu\text{m}$
- Detector: YAG crystal (0.17 mm) + Navitar zoom lens + SONY CCD camera @ 37.4 m

ADVANCED PHOTON SOURCE

A typical undulator-based emittance measurement

User run (9/22/99 4:00 AM)



Horizontal emittance = $353 \mu\text{m} \times 26.8 \mu\text{rad} \times 0.92 = 8.7 \text{ nm}\cdot\text{rad}$

(Emittance deduced from 35-BM pinhole camera results is 8.17 nm·rad at the time)

ADVANCED PHOTON SOURCE

Summary for emittance measurement with undulators

- Undulator beam contains information about divergence and size of the electron beam. It offers the potential of measuring e-beam emittance without relying on knowledge of beta functions.
- Quantitative analysis shows that it is not trivial to extract emittance information from undulator measurements. Two measurements are not enough to achieve total independence from the knowledge of lattice functions (Twiss parameters).
- We have implemented simultaneous measurement of beam "divergence" and "size" in Beamline 35-ID using a thin crystal monochromator and pinhole camera.
- Strong Fresnel diffraction characteristics have been observed in 35-ID pinhole images.


Instrument resolution compared to beam properties

Quantity	Typical value	Resolution	Comment
Horizontal divergence	26 μrad	$\sim 2.6 \mu\text{rad}$	Resolution correction not needed
Horizontal source size	340 μm	$\sim 40 \mu\text{m}$	Resolution correction not needed
Vertical divergence	5 μrad	$\sim 2.6 \mu\text{rad}$	Some resolution correction
Vertical source size	20 μm	$\sim 40 \mu\text{m}$	Resolution limited


- Improvement of vertical resolution (with imaging x-ray optics) is clearly needed.
- A third screen is also desirable in deducing emittance without relying on knowledge of Twiss parameters.

ADVANCED PHOTON SOURCE

(B) Relative beam energy measurements with undulators



Take derivative of the spectrum near the resonance, and plot against $\Delta\omega/\omega$. The result is a sharp peak centered at the nominal harmonic energy. Its FWHM is $\sim 1/nN$ (n = harmonic #, N = # of periods). The following example is calculated for the APS Diagnostics Undulator (fundamental).




ADVANCED PHOTON SOURCE

Electron Beam Energy Measurement Data

1. Scan monochromator θ - 2θ , record angle-integrated spectrum.
2. Convert the change of photon energy ω to electron energy

change:
$$\frac{\Delta p}{p_0} = \frac{\Delta \omega}{2\omega_1} = -\frac{\Delta \theta}{2 \tan \theta_1}$$

3. Fit the derivative of the spectrum with the convolution of Gaussian and SINC function and obtain the energy centroid (relative) and energy spread.



Advantages of using the undulator angle-integrated spectrum

- **Clean:** independent of beam emittance and lattice functions.
- **Accurate:** only monochromator angle needs good calibration.
- **Simple** data collection and treatment.
- **Efficient:** do not lose photons in apertures, good S/N ratio.

ADVANCED PHOTON SOURCE

(C) Energy-dependent beam property measurements

Electron beam: User fill, 100 mA in $6+25 \times 3$ bunches

Experimental Approach: Scan storage ring rf frequency (electron energy), and record

- BM and ID pinhole camera image size and centroid
- ID monochromator beam image size and centroid
- ID angle-integrated photon flux with monochromator scan
- Transient images with BM gated camera on kicker excited beam

Physical quantities extracted from the data

- Electron energy as a function of rf frequency shift → **momentum compaction factor** *
- Beam centroid as a function of beam energy → **dispersion** (x & y)
- Beam direction as a function of beam energy → **dispersion slope** (x & y)
- **Horizontal emittance** as a function of energy *
- Longitudinal emittance as a function of energy: **energy spread***
- Longitudinal emittance as a function of energy: **bunch length** (planned)
- **Transverse damping*** as a function of energy (gated camera)
- **Vertical emittance** as a function of energy (divergence only)


ADVANCED PHOTON SOURCE

Momentum Compaction Factor Measurement

- At each rf frequency shift setting
 1. Scan monochromator θ - 2θ , record angle-integrated spectrum.
 2. Fit the derivative of the spectrum with a Gaussian function.
The centroid energy is the undulator ω_1 and proportional to γ^2 .
 3. Convert the change of undulator energy ω_1 to electron energy change: $\Delta\omega/\omega_1 = 2\Delta\gamma/\gamma$
- Fit the e-beam energy change against rf frequency change to a straight line and convert to momentum compaction factor:

$$\alpha = -\frac{\Delta f}{f_0} \bigg/ \frac{\Delta p}{p_0} = -2 \frac{\Delta f}{f_0} \bigg/ \frac{\Delta\omega}{\omega_1} = \frac{2 \tan \theta_0}{f_0} \frac{\Delta f}{\Delta\theta}$$

**APS-SR Electron Energy Change with RF Frequency
(35-ID Monoscan, 5/29/99)**




For measurements of momentum compaction factors, the undulator data is very clean. One only needs to calibrate f and θ .

ADVANCED PHOTON SOURCE

Electron Beam Energy Spread Measurement

- At each rf frequency shift setting
 1. Scan monochromator θ - 2θ , record angle-integrated spectrum.
 2. Convert the change of photon energy ω to electron energy change: $\Delta\omega/\omega_1 = 2\Delta\gamma/\gamma_0$
 3. Fit the derivative of the spectrum with the convolution of Gaussian and SINC function and obtain the energy spread.




Resolution of this technique is very good: Changes of energy spread under 10^{-4} can be clearly resolved!

ADVANCED PHOTON SOURCE


Simultaneously Measured Divergence and Source Size

35-ID Emittance Measurement (Horizontal) (5/29/99)



ADVANCED PHOTON SOURCE

Comparison of ID and BM measurements






It is important to apply the correct energy spread to BM source size to deduce emittance.

$$\varepsilon_x = \left\{ \sigma_x^2 - \eta_x^2 \sigma_\delta^2 \right\} / \beta_x$$

ADVANCED PHOTON SOURCE

Transverse Damping Time (gated camera measurements)



The damping time measurements directly confirms the manipulation of the damping partition.

The shortened damping time (5.5 ms @ 800 Hz, 30% lower from 0 Hz) would further improve the performance of top-up operation if injection transient is present.

ADVANCED PHOTON SOURCE

(D) Observation of the lowest effective emittance in APS ID @ 7 GeV

Two methods to reduce emittance

- **Introduce dispersion in straight section** (reduce quantum excitation, Ed Crosbie APS/IN/ACCPHY/95-1, and Louis Emery, 5/99).

Dispersion contributes to beam size and divergence in straight section

$$\sigma_x^2 = \beta_x \epsilon_x + \left(\eta_x \frac{\sigma_E}{E} \right)^2, \text{ and } \sigma_{x'}^2 = \gamma_x \epsilon_x + \left(\eta'_x \frac{\sigma_E}{E} \right)^2$$

- **Manipulate damping partition to increase horizontal damping**

The longitudinal and transverse damping rates obey the Robinson sum rules:

$$\frac{4}{\tau_\gamma} = \frac{1}{\tau_x} + \frac{1}{\tau_y} + \frac{1}{\tau_s} = \frac{1}{\tau_\gamma} (J_x + J_y + J_s)$$

Changing partition ratio from 1:1:2 to 2:1:1 reduces horizontal emittance by half.

- Changing rf frequency shifts the electron orbit in quadrupoles in the high dispersive sections, thus changing the partition number
- Moving quadrupoles horizontally also changes the partition numbers

Applying rf frequency shift to the low-emittance lattice, we observed the lowest effective emittance in the APS storage ring to date.



ADVANCED PHOTON SOURCE

Observation of the lowest effective emittance in APS ID @ 7 GeV

Machine study run (5/29/99)

Dispersion in straight section (Crosbie-Emery low emittance lattice)

RF frequency shift +800 Hz (~3 mm orbit shift in high dispersion straight)



$$\text{Effective horizontal emittance} = 223 \mu\text{m} \times 17.1 \mu\text{rad} \times 0.92 = 3.5 \text{ nm} \cdot \text{rad}$$

ADVANCED PHOTON SOURCE

Summary for undulator based e-beam measurements

- The development of the diagnostics undulator beamline has enabled us to measure the beam emittance somewhat independent of beta function.
- Many systematic errors of the undulator-based emittance measurement have been studied theoretically and used to correct the experimental data.
- We developed a new technique for extracting e-beam energy information from undulator spectrum. It has shown good promise to become a high-resolution tool.
- The experiment for energy-dependent e-beam property measurement has illustrated the power of the undulator as a diagnostics device.
- The diagnostics undulator beamline has been instrumental in our exploration of a new operating regime of low beam emittance. One can comfortably use the data taken at this beamline to predict the performance of user ID beamlines in the low-emittance regime.

Undulator beamline upgrade

- The procurement of a cryogenically cooled monochromator is in progress. It will improve the resolution of divergence measurements with the use of the undulator third harmonic, and improve the temporal resolution of the measurement with the increased photon flux.

ADVANCED PHOTON SOURCE

IV. Summary and Future Development

- Our basic photon monitor, the x-ray pinhole camera at 35-BM has matured with the storage ring: Its stability and resolution are compatible with the current user requirements.
- Many time-resolved imaging techniques, using streak and gated cameras, have been made available to machine studies and have become critical tools for characterizing accelerator beam dynamics.
- The undulator-based beam diagnostics have gone through a rapid development process in FY1999. We expect further development in FY2000.
- Continued vigorous research and development efforts in x ray-based diagnostics are needed for the future needs of this light source and the next generation of light sources.